

Tabinfo Kuala Lumpur in November 2005

Yesterday the idea, today the sample, tomorrow the product –under this slogan the team from Heinrich Burghart GmbH presented itself at the Tabinfo in Kuala Lumpur.

This exhibition was another opportunity for us to present our well- known devices and our new developments to you.

Three new developments from our production line prodded much interest:

Particular attention was applied to the **RM 20 B** Smoking-Machine Burghart. As early as twenty years ago, Heinrich Burghart GmbH built the smoking machine RM20CS, which was distributed via Heinrich Borgwaldt GmbH. Because of many customer requirements we have revived the old RM20CS with our new **RM 20 B**.

This new Smoking-machine has been developed with new software and a touch display to satisfy the demands of the cigarette industry.

Following methods can be selected directly using the touch display:

ISO 3308, the method of Massachusetts, by the Canadian method as well as the Hoffman analysis.

But many other methods can also be created and recorded via the touch display.

Therefore all analysis can be reproduced at any time.

The **MSB-01** Mimic Smoker Burghart simulates human exposure to smoke by means of living cell cultures.

The **MUB-01** Masking Unit Burghart enabled a neutral review for smoke panel testing via covering the cigarette logo. In addition it is possible to cover up to 100% of the label and the ventilation zone.

The Tabinfo in Kuala Lumpur include a range of impressions, interesting conversations and contacts. So after three days at the exhibition we completed our way back home full of expectations and freshly motivated.